

DOLTRON-PIM 717





SERVICE MANUAL



DOLTRON AG

Industriestr. 3, CH-8610 Uster
Telefon 01 940 60 59/58
Telex 828 553 dol-ch

SPECIFICATIONS : DOLTRON PIM 717

Flow Range : 1ml/h - 69ml/min
Volume Display : 0,1 - 999,0ml
Accuracy, Pump : $\pm 2\%$
Accuracy, IV Set : $\pm 5\%$ by a random selection
of standard IV sets
Alarm Systems : audible, visual, nurse-call
Alarm Conditions : end, air (all fluids), battery
operating failure, occlusion
Power Requirements : 220V/50Hz, 110V/60Hz
Protection : Class 1
Mark of Approval : SEV, IEC 601 Part 1,    
VDE 0750 T.1/05.82
Battery : 2x12V, 1,8Ah sealed lead acid
Battery Operation : Typical operating time for fully
charge new batteries;
30ml/h ca. 3 hours
60ml/min ca. 1 hour
Case Dimensions : 170x255x220mm (WxHxD)
Weight : 5,2kg without batteries
6,9kg with batteries



DOLTRON AG TEL. 01.940 60 59-58
INDUSTRIESTR. 3 TLX. 828 533 DOL-CH
CH-8610 USTER FAX. 01.940 62 54
SCHWEIZ SWITZERLAND

15.12.1981

08.08.1985

OPERATING INSTRUCTION

Connect pump to a grounded socket outlet, unless battery operation is desired. Red dot on the panel display indicates connection to mains.

Preparing the I. V.

Consult I. V. set instruction. *Clear air completely.* Place infusion set in the pump and close lid firmly. Set prescribed ml dose-limit-total on Thumb-wheel Selector. Select desired infusion rate ml/h or ml/min on delivery rate Thumbwheel Selector. (Blinker signals ml/min operation). Make venipuncture.

Infusion

Push Green ON/OFF Button to turn on power. Push Yellow START Button to start pumping. Infused volume is on running display. Pump stops when total dose volume has been pumped.

Changing or Discontinuing I. V.

For intermediate stops, push Red STOP Button. Infused volume reading holds on display. If desired adjust Volume and/or Rate Selectors to new values. Changing of the Rate during running infusion automatically stops the pump.

CAUTION: Do not open pump lid unless clamp has been closed.

Battery Operation

The Green ON/OFF Button remains dark during battery operation.

Battery Recharge

With the pump connected to the mains supply power, the batteries will recharge regardless of operation.

Maintenance

Use cotton sticks dampened with alcohol to keep the two tubeholding sensors underneath the pump lid clean at all times. Unclean sensors will activate the Alarm (8) and stop the pump. Disconnect pump from mains power before cleaning housing. Use cloth dampened in soapy water.

ALARM INDICATIONS

By Alarm Indication the pump stops immediately. The Red STOP Button lights up and an audible alarm signals the attendant. The alarm signal may be suppressed temporarily for approx. one minute by pushing button on the back of the pump. With the remote Nurse Call Alarm connected (socket on the back of the pump), no other audible alarm will sound.

Dose Limit

The pump stops when preselected total volume and infused volume on display tally.

Air-In-Line (8)

An (8) will flash alternately with the infused volume on display. The two tube holders on each side of the peristaltic under the pump lid are both equipped with infra-red sensors for air detection in all liquids regardless of optical density.

Operational Error

Opening of the pump lid during infusion will trip the alarm.

Pump Malfunction

By malfunction, the pump stops and cannot be restarted.

Battery Alarm (b)

Estimated Battery Operation Capacity is approximately 3 hours at 30 ml/h or one hour at 60 ml/min.

With increasing battery discharge the alarm is tripped. A (b) will flash alternately with infused volume on display. Pump will continue pumping, but must be connected to mains supply power with a minimum of delay.

CHECKLIST : DOLTRON PIM 717

Alarm and Control Functions

1. Pump Lid Alarm: Start pump and open lid. Pump must stop and activate audible alarm.
2. Air-in-line : Air-Sensor control test. Start pump. Interrupt infrared light beam in the left tube holding Air-Sensor under the pump lid. Repeat for the right sensor. Pump must stop and activate audible and visual «8» alarms.

Air-Sensor control test with empty tubing.
Insert fresh, transparent piece of empty PVC-tubing (3x4mm) into the pump and start. Pump must stop and activate audible and visual «8» alarms.

Air-Sensor control test with filled tubing.
The Air-Sensors will detect air bubbles entering the pump regardless of liquid optical density. Air bubbles should trip the audible and visual «8» alarms throughout the normal infusion rate flow range viz. up to approx. 15ml/min. However, the Air-Sensors are subject to an electronically controlled time delay in order to pass air bubbles of approx. 5-8mm in the upper transfusion rates. Check at 69ml/min.
3. Dose Limit : Select a three digit dose limit between 100-999ml. Start pump. Pump must stop and activate the audible alarm when dose limit and display tally.
4. Low Battery : When the batteries run low, the audible and visual alarms are activated and signal recharging «b». The pump will continue normal operation until battery power is exhausted. Check battery power with DVM connected to L1 and L3 on board no. 25200. Battery alarm must be activated when voltage is lower than 22.0V to 22.8V. Trim with R7. Pump must cease to operate, when voltage is lower than 19.0V to 20.5V. A complete recharge of low batteries requires 14 hours.
5. Remote Alarm : The pump can be connected to the hospital nurse-call system for remote monitoring of all alarm conditions. Use a Lemo plug no. F0S304NB042 for connection to pump. Consult wiring circuitry on the back of the pump or on print no. 25200. (Doltron part no.: 25012).
6. Change in Rate: Radical changes in the delivery rate should stop the pump. No alarm.
7. Motor Control : The motor is electronically controlled in order to prevent a run away, over flow in case of a defect. Check by disconnecting Motor Control tachometer. Start pump. Pump should stop after max. 7 motor revolutions, not to be started again. 2

CHECKLIST : DOLTRON PIM 717

Alarm and Control Functions contd.

8. Occlusion : The occlusion alarm function is expressly designed to detect undue pressure build-up in the infusion range of the PIM 717 viz. 1ml/h - 30ml/min. With increasing rate sensitivity is reduced electronically in order to avoid false alarms in the transfusion range.

Check by clamping tube approx. 30 cm downstream from pump. Depending on the length of the tube enough pressure should build up to trip the alarm inside 2 - 5 ml. Pump must stop and activate audible and visual «O» alarms.

SPECIAL CHARACTERISTICS

- ON/OFF Button : If turned off, the pump can only be restarted after a delay of approx. 2 - 3 seconds.
- During battery operation, the ON/OFF Button remains dark.
- START Button : After the pump has been turned on, the electronic systems need time (max. 3 seconds) in order to coordinate prior to operation. If the START Button is pressed prematurely, the pump returns to STOP.
- Time Counter Holder : On board no.: 25200 the PIM 717 has a holder (similar to a fuse holder) intended for a cumulative electrolysis operating time counter. Not supplied. (Doltron part no.: 25061).

15.12.1981

SERVICE PROCEDURE : DOLTRON PIM 717

Prior to any service adjustments, run pump for about 10 minutes in order to warm it up.

A) Air-Sensors : print no.: 25210

1. Closed pump lid. No tubing inserted. Beware of influence from surrounding light.
2. Short M9 to Ground.
3. Connect M7 (-) and M8 (+) to DVM.
4. Balance with P4 to 0V.
5. Interrupt light beam, right sensor.
Voltage should drop to -13,2V to -14,0V.
6. Interrupt light beam, left sensor.
Voltage should rise to +13,2V to 14,0V.
7. Insert fresh, transparent piece of empty PVC-tubing.
Balance voltage with P4 to + 1.6 - 2.0V
8. Remove tubing. Voltage should drop to 0,1V to 0.5V

CAUTION: The light-sensors must be replaced in pairs.
Do not replace single sensors. Use original
by Doltron calibrated sensor pairs only. Part no. 25060
Consult color code on p. 15.

B) DA/AD converter : print no.: 25210

METRIC VERSION

1. Press STOP Button.
2. Set rate on thumbwheel to 00ml/(min or h).
3. Open bridge B2. Connect DVM at M5.
4. Set to 0,000V \pm 0,001V with P1.
5. Shorten B2 and set to 0,000V \pm 0,003V with P2.
6. Set rate on thumbwheel to 69ml and check voltage at M5.
Voltage should rise to 3,38V. For tolerances, see p. 5.
7. Connect Period-Counter at M6. Set period with P3 to 169,8 micro seconds.
8. Set rate on thumbwheel to 01ml and period duration with P2 to 11.500 micro seconds.
9. Repeat items 7 and 8!

Inch Version on
p. 4.2

Re 82.1


INCH VERSION:

Control print no. 25210

PIM 717

- 1) To be used in pumps
with the Type Nos.:

25005	110V	INCH
25006	220V	INCH
25007	240V	INCH
- 2) On board no. 25210 the PIM 717 inch-version has an "Inch-Adapter" no. 25232, instead of metric version Ic 26. The value of R57=10K



C) DA/AD converter : print no.: 25210

1. Press STOP Button.
2. Set rate on thumbwheel to 00ml/min or h.
3. Open bridge B2. Connect DVM at M5.
4. Set to $0,000V \pm 0,001V$ with P1.
5. Shorten B2 and set to $0,000V \pm 0,003V$ with P2.
6. Set rate on thumbwheel to 69ml and check voltage at M5. Voltage should rise to 3.38V. For tolerances, see p.5.
7. Connect Period-Counter at M6. Set period with P3 to 113.2 micro seconds.
8. Set rate on thumbwheel to 01ml and period duration with P2 to 7.800 micro seconds.
9. Repeat items 7 and 8 !


Re 85.0

Tolerance: $\pm 2,5\%$


Reference voltage: 15.00V (M5, print 25200 B)

Rate ml/min	Voltage at M5 (mV)			Period duration at M6 (micro sec.)		
	Nom.	Min.	Max.	Nom.	Min.	Max.
1	49	47,7	50,2	11720	11427	12013
2	98	95	100	5860	5713	6006
3	147	143	150	3906	3808	4003
4	196	191	200	2930	2856	3003
5	245	238	251	2343	2284	2401
6	294	286	301	1953	1904	2001
7	343	334	351	1674	1632	1715
8	392	382	401	1464	1427	1500
9	441	430	452	1302	1269	1335
10	490	478	502	1172	1142	1201
20	980	955	1004	586	571	600
30	1470	1433	1506	391	381	400
40	1960	1911	2009	293	286	300
50	2450	2389	2511	234	228	240
60	2940	2867	3013	195,3	190,4	200
69	3381	3296	3465	169,8	165,5	174

Layout power supply print 25200 B
BESTUECKUNGSZEICHNUNG STROMVERSORGUNGS
DOLTRON 25200 B ED 068-81.1



Layout control print 25210 B
 BESTUECKUNGSZEICHNUNG STEUERPRINT PIM 111
 DOLTRON 25210 B ED 069-81




Re 82.1

Layout display print 25220 B

BESTUECKUNGSZEICHNUNG ANZEIGEPRINT PIM 117


DOLTRON 25 220B ED 010-81




DOLIRON PIM 717 Re 10.1.82 / 20.11.82
 Bestückungszeichnung
 Layout print 25215A/25226A/ 25225/ 25227

Re 81/82.2


Muting Print 25215A




Blink Print 25226A



Blink Print 25225




Summer Print 25227




Prinzipschema PIM 717

Re April 82



Doltron PIM

Stromversorgungsprint 25200 b Re Feb.81 / Juni 81.1
Power supply print 25200B



Parts list PIM 717 power supply print 25200 B

C 1	3300mF/40V	R 19	10K	1%
C 3	100mF/35V	R 20	10K	1%
C 4	100mF/35V	R 21	150E	1%
C 5	100mF/35V	R 22	100K	1%
C 6	0,1mF/63V	R 23	56K	1%
C 7	0,1mF/63V	R 24	150E	1%
C 8	0,1mF/63V	R 25	100K	1%
		R 26	27K	1%
D 0	BZX 79C4,7	R 28	2,2M	5%
D 1	BYX 72	R 29	4,7K	1%
D 2	1N 4448	R 30*	0E Summer adjustment	
D 3	1N 4448	R 31	4,7K	1%
D 4	1N 4007	T 1	BD 234C	
D 5	1N 4007	T 2	BC 547B	
D 6	BZX 79C18	T 3	BC 547B	
D 7	1N 4448	T 4	BC 547B	
D 8	1N 4448	T 5	BC 547B	
D 9	1N 4448			
Ic 1	LM 317T			
Ic 2	LM 340T15			
Ic 3	TL 072			
RD 1	LM 336/2,5V			
R 1	10K	1%		
R 2	470E	1%		
R 3	10E	1%		
R 4	0,1E / 1W			
R 5	470E / 3W			
R 6	470E / 3W			
R 7	4,7K - 5,6K	1%		
R 8	12K	1%		
R 9	10K	1%		
R 10	10K	1%		
R 11	1,8K	1%		
R 12	1M	1%		
R 13	27K	1%		
R 14	4,7K	1%		
R 15	10K	1%		
R 16	56K	1%		
R 17	100K	1%		
R 18	1K	1%		

* Not applicable with
D 0 available

Re 82.0

Teilschema A Steuerprint 25210 b Re Feb. 81/Juni 81/Nov. 81
Part schematic A, control print 25210 B



Doltron PIM

Teilschema B Steuerprint 25210 b Re Feb.81/Juni 81




Part schematic C, control print 25210 B

Doltron PIM


Teilschema C Steuerprint 25210b

Re Feb. 81/ Juni 81/ Dez. 81/ Nov. 82




Ic 34
T 21
R 139-142
C 46-48
D 57-61

Muting-Print 25215




Ic 38
T 22
R 170-174
C 60-63

buzzer -Print 25227




Ic 40-41
T -
R 179-180
C 66
D 77

Blink-Print 25225



R 175
C 64
D 75

Blink-Print 25226



Ic 39
R 176-178
C 65
D 76

Parts list PIM 717 control print 25210 B

C 10	0,01/63V	D 10	1N 4004
C 11	0,1/63V	D 11	BZY 79C11
C 12	0,1/63V	D 12	BZY 79C11
C 13	0,01/63V	D 13	BZY 79C11
C 14	0,1/63V	D 14-44	1N 4448
C 15	0,1/63V	D 45-46	BZY 79C2,4
C 16	22mF/25V Tantal	D 47-61	1N 4448
C 17	4,7nF/63V		
C 18	1mF/25V Tantal	Ic 10	DAC 0800
C 19	0,22mF/63V	Ic 11	DAC 0800
C 20	0,22mF/63V	Ic 12	TL 072
C 21	10mF/25V Tantal	Ic 13	TL 072
C 22	2,2mF/25V Tantal	Ic 14	TL 072
C 23	0,22mF/63V	Ic 15	LM 331
C 24	0,22mF/63V	Ic 16	74 C90
C 25	0,047mF/63V	Ic 17	74 C90
C 26	2,2nF/63V	Ic 18	CD 4040
C 27	0,1mF/63V	Ic 19	74 C00
C 28	0,22mF/63V	Ic 20	74 C00
C 29	0,1mF/63V	Ic 21	74 C00
C 30	2,2nF/63V	Ic 22	CD 4098
C 31	2,2mF/25V	Ic 23	74 C193
C 32	0,1mF/63V	Ic 24	74 C08
C 33	10mF/63V	Ic 25	74 C08
C 34	22mF/25V	Ic 26	74 C90
C 35	22mF/25V	Ic 27	74 C90
C 36	0,1mF/63V	Ic 28	74 C90
C 37	0,22mF/63V	Ic 29	74 C90
C 38	10mF/25V	Ic 30	74 C90
C 39	2,2nF/63V	Ic 31	ICM 7555
C 40	2,2nF/63V	Ic 32	TL 072
C 41	0,01mF/63V	Ic 33	LM 324
C 42	2,2mF/25V	Ic 34	CD 4098
C 43	10mF/25V		
C 44	10mF/25V	T 10	BC 547B
C 45	0,01mF/63V	T 11	BC 547B
C 46	0,1mF/63V	T 12	BC 547B
C 47	2,2mF/25V	T 13	BC 547B
C 48	0,22mF/63V	T 14	BC 547B
		T 15	BC 547B
		T 16	BC 547B
		T 17	BC 547B
		T 18	BDX 33A
		T 19	BDX 33A
		T 20	BC 557B
		T 21	BDX 33A

Re 82.0

Parts list PIM 717 control print 25210 B

R 40	10K	1%		R 91	10K	
R 41	10K	1%		R 92	2,2E 4W	
R 42	10K	1%		R 93	47K	1%
R 43	5,6K	1%		R 94	470K	1%
R 44	10K	1%		R 95	10K	1%
R 45	10K	1%		R 96	47K	1%
R 46	5,6K	1%		R 97		
R 47	100K	1%		bis		
R 48	10K	1%		R 104	10K	1%
R 49	100K	1%		R 105	1M	1%
R 50	10K	1%		R 106	10K	1%
R 51	100K	1%		R 107		
R 52	100K	1%		bis		
R 53	100K	1%		R 113	100K	1%
R 54	10K	1%		R 114	22K	1%
R 55	4,7K	1%		R 115	1K	1%
R 56	100K	1%		R 116	10K	1%
R 57	5,6K	1%		R 117	10K	1%
R 58	10K	1%		R 118	10K	1%
R 59	10K	1%		R 119	10K	1%
R 60	5,6K	1%		R 120	220K	1%
R 61	10K	1%		R 121	100K	1%
R 62	10K	1%		R 122	100K	1%
R 63	10K	1%		R 123	10K	1%
R 64	10K	1%		R 124	150K	1%
R 65	10K	1%		R 125	100K	1%
R 66	10K	1%		R 126	47K	1%
R 67	10K	1%		R 127	47K	1%
R 68	10K	1%		R 128	10K	1%
R 69	10K	1%		R 129	1K	1%
R 70	100E	1%		R 130	100K	1%
R 71	100E	1%		R 131	10K	1%
R 72	10K	1%		R 132	47K	1%
R 73	10K	1%		R 133	10K	1%
R 74	10K	1%		R 134	100K	1%
R 75	2,2K	1%		R 135	100K	1%
R 76	10K	1%		R 136	100K	1%
R 77				R 137	10K	1%
bis				R 138	10K	1%
R 85	10K	1%		R 139	100K	1%
R 86	820E	1%	(1,5K)	R 140	10K	1%
R 87	820E	1%	(1,5K)	R 141	10K	1%
R 88	100K	1%		R 142	100K	1%
R 89	1K	1%		P 1	50K	
R 90	4,7K	1%		P 2	10K (new 2k)	
Ra	22K	1%		P 3	5K	
Rb	18K	1%		P 4	5K for black sensor	
				P 4	50K for brown sensor	
				P 4	100K for red sensor	
				P 4	20K for "DOLTRON Sensor"	
Re 82,2						

Connector

Stecker

K

I

H

G

F


E

D

C

B

A



Doltron PIM

Display print 25220 B
Anzeigeprint 25220b

Re Feb. 81/ Juni 81

Ic 35-37 D65-66 R 150-162 C 50

Parts list PIM 717 display print 25220 B

C 50 10mF/25V

D 65 1N 4448

D 66 1N 4448

Ic 35 74 C85

Ic 36 74 C85

Ic 37 74 C85

R 150

bis

R 162 10K 1%

Parts list PIM 717 muting print 25215

C 60 0,1mF/63V

C 61 10mF/25V

C 62 0,1mF/63V

C 63 22mF/25V

Ic 38 ICM 7555

R 170 5,6K

R 171 10K

R 172 10K

R 173 2,2M

R 174 10K

T 22 BC 547B

Parts list PIM 717 blink print 25225

C 64 100mF/25V

D 75 LED FRL-4403

R 175 1,5K

Parts list PIM 717 blink print 25226

C 65 10mF/25V

D 76 LED

Ic 39 ICM 7555

R 176 1,5K

R 177 22K

R 178 820E

Re 22.12.81 .1

Parts list PIM 717 buzzer print 25227

C 66 10mF/16V

D 77 1N 4448

Ic 40 ICM 7555


Ic 41 78L 12A

R 179 1,5K

R 180 22K

B 3 change-over switch

S buzzer, Star CMB 12



Doltron PIM 717
 Occlusion - Print 25230
 Re April 83

$U_{M1} \ 2.20V \approx t_{M2} \ 2.70ms$


Parts list PIM 717 occlusion print 25230 C

C1	10mF/25V	R11	10K/1%
C2	4,7mF/25V	R12	10K/1%
C3	1mF/25V	R13	10K/1%
C4	10nF/1%	R14	10K/1%
C5	1mF/25V	R15	10M/5%
C6	0,1mF/63V	R16	100K/1%
C7		R17	10K/1%
bis		R18	1M/1%
C10	2,2nF/63V	R19	100K/1%
C11	4,7mF/25V	R20	100K/1%
C12	22mF/16V	R21	100K/1%
C13	0,1mF/63V	R22	2,2K/1%
		R23	2,2K/1%
D1	1N4004	R24	2,2K/1%
D2			
bis		Ra	(820K/1%)
D5	1N4448	Rb	(680K/1%)
D6	LED rot	Rt	(8,2K/1%)
D7	LED gelb		
D8	LED grün		
Ic A	LM331N		
Ic B	CD4040		
Ic C	CD4040		
Ic D	CD4011		
Ic E	CD4011		
Ic F	CD4040		
Ic G	CD4011		
Ic H	ICM7555		
Ic I	CD4585		
Ic K	CD4585		
Ic L	CD4040		
Ic M	CD4011		
R1	100K/1%		
R2	470K/1%		
R3	51K/1%		
R4	10K/1%		
R5	27K/1%		
R6	10K/1%		
R7	100K/1%		
R8	1K/1%		
R9	10K/1%		
R10	10K/1%		

Re April 83.1

Layout occlusion print


25230 C




DOLTRON Re April 83

PIM 717 Steuerprint 25210 Lötseite

Connections for occlusion print 25230 C




Re April 83.1



Parts list PIM 717
Soft-Start Print 25231C

C67	1 mF/16V
D68	Ref. D. 6,2V
R180	8,2E/4W
R181	3,9E/4W
R182	100K/1%
R183	10K/1%
R184	270K/1%
R185	10K/1%
R186	4,7M/5%
R187	10K/1%
R188	470E/1%
REL.3	24V/1xU
Ic42	LM 358N



Doltron PIM 717


Soft-Start Print 25231C

Re Juni 83

Type 305

7-SEGMENT DISPLAY WITH BCD INPUT AND MEMORY

Circuit diagram



In this module, the display is operated by the BCD code in positive logic. A control signal (Input EL) makes it possible to freeze the display and suppress response to changing BCD input signals.

Caution: The terminals for TTL and C-MOS are not identical.

Technical data

Colour of display	bright red
Digit height	9 mm
Supply current I_{cc} for TTL	typ. 93 mA
for C-MOS	" 48 mA

Input data

Input voltage (all inputs)	TTL	C-MOS
U_{in} "0" with $V_{cc} = 5V$	max. 0.8 V	1.5 V
$= 10V$	max.	3 V
$= 12V$	max.	3.6 V
$= 15V$	max.	4.5 V
U_{in} "1" with $V_{cc} = 5V$	min. 2 V	3.5 V
$= 10V$	min.	7 V
$= 12V$	min.	8.4 V
$= 15V$	min.	10.5 V

Input currents

TTL logic: BCD input currents with EL = "0"

I_{in} "0" with $U_{in} = 0.4 V$	max. - 1.6 mA
I_{in} "1" with $U_{in} = 2.4 V$	max. 80 uA
With EL = "1" $I_{in} =$ "0" and "1"	max. - 0.1 mA

Inputs EL and RBI

I_{in} "0" with $U_{in} = 0.4 V$	max. - 1.6 mA
I_{in} "1" with $U_{in} = 2.4 V$	max. 40 uA

Input BI

I_{in} "0" with $U_{in} = 0.4 V$	max. - 3.2 mA
I_{in} "1" with $U_{in} = 2.4 V$	max. - 80 uA

C-MOS logic: all inputs

I_{in} "0" and I_{in} "1"	typ. 10 pA
Input capacitance C_{in}	typ. 5 pF

Output data (RBO only)

Output voltage

U_{out} "0" with $I_{out} = - 3.2 mA$	max. 0.4 V
U_{out} "1" with $I_{out} = - 80 uA$	min. 2.4 V

Output current

I_{out} "0"	max. - 3.2 mA
I_{out} "1"	max. - 80 uA

Description of Enable Latch (EL) and Decimal Point (DP):

EL (Enable Latch) This instruction will freeze the display and suppress further response to changes of the BCD input.

"EL" on "0"	The display responds to the BCD input value
"EL" on "1"	The display freezes on the last value

DP (Decimal Point) The decimal point must be controlled externally. The module features an integral current limiting resistor.

"DP" on "0"	Decimal point off
"DP" on "1"	Decimal point on

Truth tables

C-MOS logic

Inputs							Output
EL	LT	D ₂ ³	C ₂ ²	B ₂ ¹	A ₂ ⁰	BI	Display
X	0	X	X	X	X	X	B (Test)
X	1	X	X	X	X	0	none
0	1	0	0	0	0	1	0
0	1	0	0	0	1	1	1
0	1	0	0	1	0	1	2
0	1	0	0	1	1	1	3
0	1	0	1	0	0	1	4
0	1	0	1	0	1	1	5
0	1	0	1	1	0	1	6
0	1	0	1	1	1	1	7
0	1	1	0	0	0	1	8
0	1	1	0	0	1	1	9
1	1	X	X	X	X	1	stored*

X = "0" or "1"

* controlled by applied BCD code during the leading edge of the "EL" instruction signal

TTL logic

Inputs							Outputs	
EL	RBI	D ₂ ³	C ₂ ²	B ₂ ¹	A ₂ ⁰	BI ¹⁾	RBO	Display
X	X	X	X	X	X	0	0	none
0	0	0	0	0	0	0	0	none
0	1	0	0	0	0	0	1	0
0	X	0	0	0	1	1	1	1
0	X	0	0	1	0	1	1	2
0	X	0	0	1	1	1	1	3
0	X	0	1	0	0	1	1	4
0	X	0	1	0	1	1	1	5
0	X	0	1	1	0	1	1	6
0	X	0	1	1	1	1	1	7
0	X	1	0	0	0	1	1	8
0	X	1	0	0	1	1	1	9
0	X	1	0	1	0	1	1	A
0	X	1	0	1	1	1	1	b
0	X	1	1	0	0	1	1	c
0	X	1	1	0	1	1	1	d
0	X	1	1	1	0	1	1	e
0	X	1	1	1	1	1	1	F
1	X	X	X	X	X	1	1	stored*


X = "0" or "1"

* controlled by the applied BCD code during the leading edge of the "EL" instruction signal

1) Input BI should only be shifted to "0" to obtain blanking of the display irrespective of the BCD input. Further information on this input is provided in the general data section.

Ordering information for series 300

Display modules:



Example: Type 04 for 10 V, black housing
Order code 304-054-21

Accessories:

- 300-01-101 = grey dummy housing
- 300-01-102 = black dummy housing
- 300-03-301 = grey end bracket pair
- 300-03-302 = black end bracket pair
- 300-04-311 = thumbwheel switch adapter, grey, left-hand mounting
- 300-04-312 = thumbwheel switch adapter, grey, right-hand mounting
- 300-04-321 = thumbwheel switch adapter, black, left-hand mounting
- 300-04-322 = thumbwheel switch adapter, black, right-hand mounting
- 300-06-.L. = contrast filter (L = length in mm 016, 024, 032)
- 300-08-.L. = mounting clip (L = length in mm)*
- 000-20-.L. = Normally not needed. For larger blockassemblies are mounting clips available.


SPAREPARTS LIST: 985, PIM 717

Art No	Designation
24050	Label, "SEV, IEC, VDE"
25031	Handle
25034	Mains power socket connector
25036	Compensating plug
25038	Micro switch
25042	Buzzer 12V
25058	Motor sensor
25060	Air sensor (pair, color coded)
25111	ON/OFF switch, excl. lamp/cap
25112	Switch START/STOP excl. lamp/cap
25113	Cap ON/OFF
25114	Cap EIN/AUS
25115	Cap START
25116	Cap STOP
25118	Lamp bulb 14V/80mA
25119	Fuse 630mA/5x20mm
25121	Cap remover (tool)
25122	Lamp remover (tool)
25126	Endpiece, pair, Multiswitch/Displ
25127	Thumbwheel switch segment
25131	Spacer "m1"
25132	Spacer "TOTAL"
25281	Bottle pole holder
25282	Clamp screw
25283	Plastic tip for 25282
25287	Lock, snap catch
25289	Spacer for no. 25034
25290	Rubber feet
25010	Lead acid batteries 2x12V/1.8Ah
25029	Pump cover
25032	Label, pump door
25043	Label, Operating Instruction ger
25044	Label, Operating Instruction fr.
25045	Label, Operating Instruction eng
25048	Transformer 110V/60Hz
25049	Transformer 240V/50Hz
25050	Transformer 220V/50Hz
25120	Fuse 2A/5x20mm
25128	Thumbwheel switch, 0-6
25129	Thumbwheel switch, m/h
25130	Spacer, "/" with LED
25133	Spacer, blank, left of display
25134	Spacer, blank, right of display
25156	LED display segment
25157	LED red filter
25270	Frontplate
25272	Backplate
25200	Power supply print
25210	Control print
25210i	Control print with inch adapter
25232	Inch adapter for no. 25210
25215	Muting print
25220	Display print
25230	Occlusion print
25215	Muting print
25226	Blink print
25227	Summer print
25231	Soft-Start print

Betr. Gerät(e): PIM 303 und PIM 717


Problem: Nicht erkannte Luftförderung bei falsch eingelegtem Infusionsschlauch.

Ursache: Der Infusionsschlauch wird nicht, wie in der Bediennungsanleitung beschrieben und dokumentiert, in einer Linie durch die Luftsensoren und die Peristaltik geführt. Es besteht die Möglichkeit, mit einigem Geschick den Infusionsschlauch um eine oder beide Luftsensoren herum zu führen. Wenn der Sensor auf der Eingangsseite (links) der Peristaltik umfahren wird, kann das durch die Elektronik nicht erkannt werden.



Lösung: 1. Klare Instruktion der Personen, die mit PIM 303 und PIM 717 arbeiten. Eventuell fehlende Bedienungsanleitungen nachbestellen.

2. Einbau von zusätzlichen Schlauchführungen im Bereich der beiden Luftsensoren. Durch diese Schlauchführungen wird das richtige Einlegen des Infusionsschlauchs verdeutlicht. Wird der Infusionsschlauch trotzdem ganz bewusst um die Schlauchführungen und Sensoren herumgeführt, so kann der Pumpendeckel nur noch mit erheblichem Kraftaufwand geschlossen werden. In diesem Fall wird der Infusionsschlauch aber abgeklemmt. Es entsteht die Situation der "Nullförderung".



Hinweis: Als vorbeugende Massnahme sollten alle Pumpen der PIM-Serie beim nächsten Service mit Schlauchführungen ausgerüstet werden. Bei Pumpen des Typs PIM 717, die zusammen mit Spezialbestecken der Serie "HIFLOW" eingesetzt werden, dürfen nur Schlauchführungen eingesetzt werden, die für den Einsatz von HIFLOW-Bestecken nachbearbeitet worden sind.

Art.Nr.:

25422	Schlauchführung "links", inkl. Senkschrauben
25423	Schlauchführung "rechts", inkl. Senkschrauben
25427	Schlauchführung "links" Typ HIFLOW, inkl. Senkschrauben
25428	Schlauchführung "rechts" Typ HIFLOW, inkl. Senkschrauben

DOLTRON Code Opto-Sensors

Used as Air Sensors

Group	Color code	CTR in %	PIM 303		PIM 717	
			Board No.: 24225		Board No.: 25210	
			P1	R1 / R2	P4	R86 / R87
1	blue/red	5-10	---	---	---	---
2	yellow/red	11-15	50K	1K	50K	1.5K
3	yellow/red	16-20	50K	1K	50K	1.5K
4	yellow/red	21-25	50K	1K	50K	1.5K
5	green/red	26-30	50K / 20K	1K	50K / 20K	1.5K
6	black/red	31-35	20K	1K	20K	1.5K
7	black/red	36-40	20K	1K	20K	1.5K
8	black/red	41-45	20K	1K	20K	1.5K
9	white/red	46-50	20K	2.2K	20K	2.7K
10	white/red	51-55	20K	2.2K	20K	2.7K
11	white/red	56-60	20K	2.2K	20K	2.7K
12	white/red	61-65	20K	2.2K	20K	2.7K

Used as Motor-Sensors

Group	Color code	CTR in %	PIM 303	PIM 717
			Board No.: 24210	Board No.: 25210
			R32	R129
1	blue/red	5-10	820E	1K
6-8	black/red	31-45	1.5K	1K
9-12	white/red	46-65	3.3K	1K

Re: PIM 303 und PIM 717

Problem: The Opto Sensor used as PIM **air-sensors** and **motor-sensors** after having been especially selected, **color coded and grouped together in pairs of blue/red, yellow/red und green/red**, are now no longer available an the market.

Reason: The efficiency (CTR) of these components has been gradually improved over the years. This has also changed their characteristics. For use our pumps we have always had to individually check and group these together in suitable pairs be used for air sensing in PIM. With the improved efficiency of the components, we today can group sensors together coded **black/red** and **white/red** only.

Solution: Please observe, that for the proper function of the PIM air detection feature, viz. motor sensors, when using these new components coded **black/red** and **white/red some resistance values of the circuit in question have to be changed accordingly.**
Exactly which components that matter can be seen in the **Code Opto-sensors document.**

Important: When performing the "Adjustment Air-Sensors", please, consult and strictly follow the respective PIM 303 and PIM 717 Service-Manuals

Hint: *The referenced documents are:*
- Code Opto-Sensors
- Service-Manual DOLTRON PIM 303
- Service-Manual DOLTRON PIM 717

Art.No.:	Air-Sensor pair	25060 + Color code
	Motor-Sensor	25058 + Color code
	Potentiometer 20K PIM 303	98301
	Potentiometer 50K PIM 303	98302
	Potentiometer 20K PIM 717	98303
	Potentiometer 50K PIM 717	98304
	Resistor 820E/0.6W	98310
	Resistor 1.5K/0.6W	98311
	Resistor 2.2K/0.6W	98312
	Resistor 2.7K/0.6W	98313
	Resistor 3.3K/0.6W	98314

Contact:	DOLTRON AG	Tel.	++41 (0)1 940 60 58
	Industriestr. 3	FAX	++41 (0)1 940 62 54
	CH-8610 Uster	e-mail	doltron@bluewin.ch